**GEM™ CVP 2 (Calibration Valuation Product) with CO-Ox - 00025000125****Intended Use**

External Calibration Valuation Product used to complete the calibration process of the GEM Premier 3500/3000 and GEM Premier 4000 analyzers prior to use with patient samples.

**Summary and Principle**

GEM CVP for use with the GEM Premier 4000 analyzer is comprised of the following four ampule configurations:  
 GEM CVP 1 with CO-Ox: Low pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, COHb, MetHb, HbHb and high pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate. HbH and O<sub>2</sub>Hb values.  
 GEM CVP 2 with CO-Ox: High pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, IBil<sup>+</sup>, Cl<sup>-</sup>, COHb, MetHb, HbH and low pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate. HbH and O<sub>2</sub>Hb values.  
 GEM CVP 3: Low hematocrit value.  
 GEM CVP 4: Normal hematocrit value.  
 GEM CVP for use with the GEM Premier 3500/3000 analyzer is comprised of the following four ampule configurations:  
 GEM CVP 1 with CO-Ox: Low pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, and high pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate values.  
 GEM CVP 2 with CO-Ox: High pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, and low pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate values.  
 GEM CVP 3: Low hematocrit value.  
 GEM CVP 4: Normal hematocrit value.

**Composition**

GEM CVP 1 and 2 with CO-Ox ampules: Aqueous buffered bicarbonate solution containing inorganic salts and organic metabolites, stabilizer, surfactant, dye and biocides; equilibrated with precise concentrations of carbon dioxide and oxygen. Each ampule contains 1.8 mL of solution.

**GEM™ CVP 2 (Calibration Valuation Product) with CO-Ox - 00025000125****Beschreibung**

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**GEM™ CVP 2 (Calibration Valuation Product) with CO-Ox - 00025000125****Utilisation**

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GEM CVP 3 and 4 ampules: Equilibrated aqueous buffered bicarbonate solution containing inorganic salts and biocides; equilibrated with carbon dioxide and oxygen. Each ampule contains 1.8 mL of solution.  
 Hazard class: **None**  
 Risk phrases: **None**  
 Safety phrases: **None**  
 For in vitro diagnosis use.

**Storage and Stability**

GEM CVP 1 and 2 with CO-Ox: Unopened ampules are stable until the expiration date stated on the label when stored at 2-8°C or up to 8 months at room temperature (15-25°C), providing storage does not exceed the expiration date. DO NOT FREEZE. GEM CVP 3 and 4: Unopened ampules are stable until the expiration date shown on the label when stored at room temperature (15-25°C, DO NOT FREEZE.

**WARNING:** GEM CVP with CO-Ox solutions are sensitive to ambient temperature variations (gas/liquid equilibrium) and room air contamination (diffusion gradients). Ampules equilibrated at temperatures cooler than 22 +/- 1°C will demonstrate PO<sub>2</sub> values HIGHER than those stated on the label by approximately 10% per °C. Likewise, GEM CVP ampules equilibrated at temperatures warmer than 22 +/- 1°C will demonstrate PO<sub>2</sub> values LOWER than those stated by approximately one percent per °C. There is no significant temperature effect on pH and pCO<sub>2</sub> within the normal range of room temperature.

**Directions for Use GEM Premier 3500/3000 and 4000 Analyzer Users**

1. Immediately prior to use, hold the ampule by the top above the break line and shake vigorously for approximately 10 seconds.
2. Restore liquid to bottom of vial with gentle tapping.
3. Carefully snap open the ampule utilizing the ampule breaker.
4. Carefully snap open the ampule utilizing the ampule breaker.

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**ENGLISH - Insert Revision 10/2012**

5. Aspirate samples immediately after opening. Any delay in measuring may cause room air contamination and result in lower pCO<sub>2</sub> and higher PO<sub>2</sub> values than stated on the insert.
6. Refer to the appropriate Operator's Manual or Quick Operator's Guide for further information regarding product use.  
 NOTE: The CVP ampules are single use only. Do not store after opening.

**Limitations**

The values obtained from GEM CVP with CO-Ox solutions are sensitive to ambient temperature variations (gas/liquid equilibrium) and room air contamination (diffusion gradients). The value solution in each ampule immediately upon opening.

**Traceability**

Ampules were equilibrated to 22°C and measured at 37°C. Established protocols for GEM CVP are traceable to NIST materials and recognized CLSI methods where available and appropriate.

**Expected Values**

The reference values for this product were based upon multiple determinations performed on GEM Premier 3500/3000 and GEM Premier 4000 systems.

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