**GEM® CVP 4 (Calibration Valuation Product) Hematocrit - 00025000140**

ENGLISH - Insert Revision 11/2014

**Intended Use**

External Calibration Valuation Product used to complete the calibration process of the GEM Premier 3500/3000 and GEM Premier 4000 analyzers prior to use with patient samples.

**Summary and Principle**

GEM CVP for use with the GEM Premier 4000 analyzer is comprised of the following four ampule configurations:

- GEM CVP 1 with CO-Ox: Low pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, COHb, MetHb, Hhb and high pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate, Hhb and O<sub>2</sub>Hb values.
- GEM CVP 2 with CO-Ox: High pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, COHb, MetHb, Hhb and low pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate, Hhb and O<sub>2</sub>Hb values.
- GEM CVP 3: Low hematocrit value.
- GEM CVP 4: Normal hematocrit value.

GEM CVP for use with the GEM Premier 3500/3000 analyzer is comprised of the following four ampule configurations:

- GEM CVP 1 with CO-Ox: Low pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, and high pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate values.
- GEM CVP 2 with CO-Ox: High pH, pO<sub>2</sub>, Na<sup>+</sup>, K<sup>+</sup>, and low pCO<sub>2</sub>, Ca<sup>++</sup>, glucose and lactate values.
- GEM CVP 3: Low hematocrit value.
- GEM CVP 4: Normal hematocrit value.

**Composition**

GEM CVP 1 and 2 with CO-Ox ampules: Aqueous buffered bicarbonate solution containing inorganic salts and organic metabolites, stabilizer, surfactant, dye and biocides; equilibrated with precise concentrations of carbon dioxide and oxygen. Each ampule contains 1.8 mL of solution.

**GEM® CVP 4 (Calibration Valuation Product) Hematocrit - 00025000140****Beschreibung**

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ENGLISH - Insert Revision 11/2014

**3. Restore liquid to bottom of vial with gentle tapping.**

4. Carefully snap open the ampule utilizing the ampule breaker.
5. Aspirate samples immediately after opening. Any delay in measuring may cause room air contamination and result in lower pCO<sub>2</sub> and higher PO<sub>2</sub> values than stated on the insert.
6. Refer to the appropriate Operator's Manual or Quick Operator's Guide for further information regarding product use.

NOTE: The CVP ampules are single use only. Do not store after opening.

**Limitations**

The values obtained from GEM CVP with CO-Ox solutions are sensitive to ambient temperature variations (gas/liquid equilibrium) and room air contamination (diffusion gradient). Be sure the solution in each ampule immediately upon opening.

**Traceability**

Ampoules were equilibrated to 22°C and measured at 37°C. Established protocols for GEM CVP are traceable to NIST materials and recognized CLSI methods where available and appropriate.

**Expected Values**

The reference values for this product were based upon multiple determinations performed on GEM Premier 3500/3000 and GEM Premier 4000 systems.

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