

## HEMOCUE HB 201 DM Competency Test Competency Exam

This exam is to be performed upon Initial Orientation, Training, and Competency, as well as, upon Annual Competency testing. (Newly trained staff must perform this test *twice* within the first year.)

**THIS TEST MUST BE PERFORMED IN FRONT OF ANOTHER CERTIFIED OPERATOR  
(CERTIFIED OPERATOR = A STAFF MEMBER CURRENTLY APPROVED TO PERFORM TEST)**

<b>Printed Name:</b>	<b>License #:</b>
<b>Location:</b>	<b>Date:</b>

**\*\*\*\* I HAVE READ THE HEMOCUE HB 201 DM POLICY AND PROCEDURE BEFORE TAKING  
THIS EXAM AS INDICATED BY MY INITIALS HERE: \_\_\_\_\_.**

WRITTEN TEST QUESTIONS – ANSWER <u>TRUE</u> OR <u>FALSE</u>	ANSWER
1. Making up a medical record number (MRN) is okay if it is an emergency.	
2. If I mistrust a result, for any reason, I should repeat the test immediately.	
3. Microcuvettes are stored at 2-8 degrees centigrade until the expiration date located on the labels.	
4. Hemocue HB 201 DM testing is considered a WAIVED test according to the FDA	
5. Instrument maintenance and Quality Controls must be performed monthly prior to patient testing.	
6. If patient result is >23.5 g/dL, a qualified personnel must be confirmed by sending a patient specimen to the Clinical Laboratory.	
7. The presence of an air bubble in the microcuvette can produce erroneously low readings.	
8. Quality controls are stored at 15-30 degrees centigrade until the expiration date located on the labels.	
9. Sharing my Bar Code ID is okay if it is an emergency.	
10. Excessive squeezing of the finger can dilute the sample with tissue fluid and may give higher results.	
You must score <b>100%</b> on the written test (above) to proceed to test demonstration (below).	SCORE

### PLEASE PERFORM 2 LEVELS OF LIQUID QUALITY CONTROL

<b>Hemoglobin LQC Level 1</b>	<b>Hemoglobin LQC Level 3</b>
_____ g/dL	_____ g/dL

<b>I Performed This Test:</b>	<b>I Observed This Test Being Performed:</b>
<b>Signature:</b> _____ <b>Date:</b> _____	<b>Printed Name:</b> _____ <b>Date:</b> _____

**UPON COMPLETION: Fax copy to POCT Services @ 415.206.3451 or 415.206.3045 or bring copy to  
2M14. Unit Manager is encouraged to retain a copy in employee or unit file.**