GEM® Premier 4000 PAK

Intended Use
Sensor and reagent cartridge for use with the IL GEM Premier 4000 analyzer.
Storage Requirement
STORE at 15°-25° C
DO NOT FREEZE
Do not open cartridge protective wrapper until ready to use.

Instructions for Use
To insert a GEM Premier 4000 PAK cartridge:
1. Unpack the cartridge from the protective wrapper. Ensure that the cartridge is dry and within expiration date prior to use.
2. Remove plastic backing from cartridge by pulling tab.
3. Press the Open Door button. The door will release. Fully open the door to the left.
Note: The cartridge must be at 15-25°C prior to use.
4. Position the cartridge with the gray sampling area facing forward. Push the cartridge in until you feel resistance.
Note: The cartridge cannot be inserted all the way into the compartment by hand. A small lip of the cartridge will rest on the door.
5. Guide the system door to the right to close it and move the cartridge into the final position. The door will lock and the warm-up cycle will take approximately 40 minutes. Samples cannot be analyzed during warm-up.
Note: Do not force door closed. If resistance is met when closing the door, please remove the cartridge and restart the process.
6. Run appropriate GEM CVP (Calibration Valuation Product) levels prior to analyzing patient samples.
Refer to the GEM Premier 4000 Operator's Manual or Operator's Guide for complete instructions.

Warning
GEM 4000 A Solution
Hazard class: None
Hazard statements: None
Precautionary statements: None
Supplemental hazard information: Up to 2.95% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
GEM 4000 B Solution
Hazard class: Aquatic chronic 3, H412
Hazard statements: H412: Harmful to aquatic life with long lasting effects.
Precautionary statements: P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.
Supplemental hazard information: EUH208: Contains Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1). May produce an allergic reaction. Up to 2.15% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
GEM 4000 C Solution
Warning
Hazard class: Sensitization - Skin, Cat. 1
Hazard statements: May cause an allergic skin reaction.
Precautionary statements: Avoid breathing vapors/spray. Wear protective gloves. IF ON SKIN: Wash with plenty of water. If skin irritation or a rash occurs: Get medical advice/attention.
Supplemental hazard information: EUH208: Contains m-Phenyldiamine. May produce an allergic reaction. Up to 2.08% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
GEM 4000 D Solution
Hazard class: None
Hazard statements: None
Precautionary statements: None
Supplemental hazard information: GEM 4000 D Solution: EUH208: Contains Gentamicin sulfate (salt). May produce an allergic reaction. Up to 1.4% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

ENGLISH - Insert Revision 09/2016

GEM 4000 Reference Solution
Warning
Hazard class: Eye Irrit. 2, H319; Aquatic chronic 3, H412
Hazard statements: H319: Causes serious eye irritation. H412: Harmful to aquatic life with long lasting effects.
Precautionary statements: P280: Wear eye protection/face protection. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313: If eye irritation persists: Get medical advice/ attention. P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.
Supplemental hazard information: 2.22% of the mixture consists of ingredients of unknown acute toxicity (inhalation) for human health.
GEM 4000 Lysing Solution
Danger
Hazard class: Eye dam.1, H318; Skin Sens. 1, H317, Aquatic chronic 3, H412
Hazard statements: H318: Causes serious eye damage. H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects.
Precautionary statements: P280: Wear protective gloves/protective clothing/ eye protection/face protection. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302+P352: IF ON SKIN: Wash with plenty of water. P310: Immediately call a POISON CENTER/doctor. P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.
Supplemental hazard information: Contains: Polyethylene glycol octylphenyl ether and Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1). Up to 3.5% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.
For *in vitro* diagnostic use.
The GEM® Premier 4000 PAK is protected by one or more of the following US Patents: 5,132,345, 5,286,364, 5,540,828, 6,652,720, 6,872,297, 6,960,466, 7,022,219

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