PROCEDURE/INSTRUMENT: Syphilis Health Check Competency Exam

☐ Initial Competency ☐ Semi-annual Competency ☐ Annual Competency (CY:) ☐ Fo	llow up
This exam is to be performed after Initial Orientation & Training and prior to perform	ing
technical duties, 6 months after initial competency assessment, and annually thereaf	ter.

THIS TEST <u>MUST</u> BE PERFORMED IN FRONT OF ANOTHER CERTIFIED OPERATOR (CERTIFIED OPERATOR = A STAFF MEMBER <u>CURRENTLY APPROVED</u> TO PERFORM TEST)

I HAVE <u>READ</u> THE SYPHILIS HEALTH CHECK POLICY AND PROCEDURE AND THE CHANGES THEREIN BEFORE TAKING THIS EXAM AS INDICATED BY MY INITIALS HERE: ______.

WRITTEN TEST QUESTIONS – ANSWER TRUE OR FALSE	ANSWER
1. A plasma specimen may be used for testing.	
2. If I mistrust a result, for any reason, I should repeat the test immediately.	
3. The Syphilis Health Check test is considered a High Complexity test according to the FDA.	
4. The kit is stable until the expiry date stated on the package label.	
5. Dispense 25 uL of whole blood into the sample well.	
6. If there is no color line visible in the "C" control area, whether or not there is a line in the	
"T" test area, the test in invalid and cannot be interpreted.	
7. The patient result may be read any time before 3 minutes or after 5 minutes.	
8. Only 1 line must appear on the kit in order to read the patient result as positive.	
9. It is acceptable to use the first drop of blood for testing.	
10. If a test fails QC, it is considered invalid and must be repeated.	
CTON	SCORE
STOP!	
You must score 100% on the written test (above) to proceed to test demonstration (below).	
, , , , , , , , , , , , , , , , , , , ,	

Printed Name:	License # (Required):	Location:
Syphilis Health Check Kit Lot #:	Kit Exp. Date:	Blind Sample #:
1. QC line was present in control window:	YES	NO
QC background cleared or remained clear:	YES	NO
3. "Patient" (Blind Sample) Result:		

I Performed This Test (Signature):	I Observed This Test Being Performed (Signature):
/Date:	/Date: