1. I just finished sampling an arterial blood gas specimen and I realized that I selected the "venous" specimen type by mistake. Circle the best answer: a. Enter a note in the COMMENT field about it and notify POCT as well as the provider. b. I immediately re-test the specimen on another GEM 4000 if available and notify POCT and the provider about the error on the first instrument. c. Both a and b are correct. 2. In preparing a syringe blood gas sample, you should: Mix the sample for 30 seconds a) alternating between inverting and rolling between outstretched palms. Then I expel a few drops of the sample onto a gauze pad in order to check for b) which can lead to c) Circle one answer in each column: (a)	Printed	Employee Name:		Location:			
specimen type by mistake. Circle the best answer: a. Enter a note in the COMMENT field about it and notify POCT as well as the provider. b. I immediately re-test the specimen on another GEM 4000 if available and notify POCT and the provider about the error on the first instrument. c. Both a and b are correct. 2. In preparing a syringe blood gas sample, you should: Mix the sample for 30 seconds a) alternating between inverting and rolling between outstretched palms. Then I expel a few drops of the sample onto a gauze pad in order to check for b) which can lead to c) Circle one answer in each column: (a) (b) (c) 1. Right after collection and 2. Hemolysis damage the GEM cartridge asmpling 2. An increase in pO2 3. Longer analysis time by the instrument 3. I scanned the patient ID barcode and checked the demographics on the screen. I realized that I used the wrong patient ID barcode. I immediately notify a) and b) Afterwards, I notify c) using the "Wrong Patient ID / Wrong Specimen Type" form which I submit via email or at the POCT office mailbox. Circle one answer in each column: (a) (b) (c) 1. The provider of the (a) (b) (c) 1. The provider of the (b) (c) 1. The provider of the (a) (b) (c) 1. The provider of the (b) (c) 2. Recollect a specimen	GEM PREMIER 4000 PROBLEM SOLVING TEST						
c. Both a and b are correct. 2. In preparing a syringe blood gas sample, you should: Mix the sample for 30 seconds a)	1.	specimen type by mistake. Circle the best answer: a. Enter a note in the COMMENT field about it and notify POCT as well as the provider. b. I immediately re-test the specimen on another GEM 4000 if available and notify POCT					
(a) (b) (c) 1. Right after collection 2. Right after collection and then right before 3. Clots 3. Clots 3. Clots 4. I scanned the patient ID barcode and checked the demographics on the screen. I realized that I used the wrong patient ID barcode. I immediately notify a) and b) Afterwards, I notify c) using the "Wrong Patient ID / Wrong Specimen Type" form which I submit via email or at the POCT office mailbox. Circle one answer in each column: (a) (b) (c) 1. The provider of the correct patient	2.	c. Both a and b are correct. In preparing a syringe blood gas sample, you should: Mix the sample for 30 seconds a) alternating between inverting and rolling between					
1. Right after collection 2. Right after collection and then right before sampling 3. Clots 3. Clots 4. An increase in pO2 3. Longer analysis time by the instrument 3. I scanned the patient ID barcode and checked the demographics on the screen. I realized that I used the wrong patient ID barcode. I immediately notify a) and b) Afterwards, I notify c) using the "Wrong Patient ID / Wrong Specimen Type" form which I submit via email or at the POCT office mailbox. Circle one answer in each column: (a)		for b) which can lead to	o c)	Circle one answ	le one answer in each column:		
3. I scanned the patient ID barcode and checked the demographics on the screen. I realized that I used the wrong patient ID barcode. I immediately notify a) and b) Afterwards, I notify c) using the "Wrong Patient ID / Wrong Specimen Type" form which I submit via email or at the POCT office mailbox. Circle one answer in each column: (a)		 Right after collection Right after collection and then right before sampling 	2. Hemolys	, ,	damage t cartridge 2. An increa 3. Longer ar	esults and the GEM ase in pO2 nalysis time by	
1. The provider of the correct patient. 2. The provider of the wrong patient ID I just used. 3. The patient's next of kin. 4. I scanned the CSN barcode into the ACCOUNT field and no demographics displayed on the GEM screen. Following the procedure, I verified that I used the correct patient ID. Next, I must enter the in the PATIENT ID field. Circle the best answer: a. Patient's initials. b. Patient's last name. c. Patient's MRN number. d. Last four digits of the CSN number. YOU MUST SCORE 100% TO PASS. If you did not pass, please review the Policy and Procedure and retake the exam.	3.	used the wrong patient ID barcode. I immediately notify a) and b) Afterwards, I notify c) using the "Wrong Patient ID / Wrong Specimen Type" form					
screen. Following the procedure, I verified that I used the correct patient ID. Next, I must enter the in the PATIENT ID field. Circle the best answer: a. Patient's initials. b. Patient's last name. c. Patient's MRN number. d. Last four digits of the CSN number. YOU MUST SCORE 100% TO PASS. If you did not pass, please review the Policy and Procedure and retake the exam.		 The provider of the correct patient. The provider of the wrong patient ID I just used. 	 Recollect from the and re-te Notify th 	t a specimen correct patient est. e nurse	2. The nurse	vices e educator.	
If you did not pass, please review the Policy and Procedure and retake the exam.	screen. Following the procedure, I verified that I used the correct patient ID. Next, I must enter the in the PATIENT ID field. Circle the best answer: a. Patient's initials. b. Patient's last name. c. Patient's MRN number. d. Last four digits of the CSN number.						
Employee Signature / Date Trainer Signature / Date	100 W031 3COKE 100% TO FA33.						
	Employee Signature / Date			Trainer Signature / Date			