

Cue™ COVID-19 Competency Exam

Initial Competency Semi-annual Competency Annual Competency (CY: _____) Follow up

This exam is to be performed upon Initial Orientation & Training and Competency; twice in the first year of practice; and annually thereafter.

PRINTED Name:	License #:	Home Unit:
I have READ/REVIEWED the procedure prior to taking this test as indicated by my initials: _____		
DATE of Test:	TIME of Test:	
Result of blind sample:	Neg or Pos (circle one)	

Written Exam – Circle either T = statement is true or F = statement is false

T	F	1. It is acceptable to use any collection swab as long as it is clean.
T	F	2. Cue Test control swabs must be tested with each new shipment and new lot of test cartridges received and when problems occur with testing, prior to patient testing.
T	F	3. An invalid result means that there was merely something wrong with the Sample Wand.
T	F	4. Within about 20 minutes of inserting the swab into the test card, the result will be available.
T	F	5. The result is saved in the Cue Account profile that was selected before the test started.
T	F	6. The Cue Health App shows the result as Negative, Positive, Invalid, or Canceled.
T	F	7. Negative results should be considered confirmed and no further testing is required even with clinical signs and symptoms.
T	F	8. You must insert the Sample Wand with nasal sample into the Cue COVID-19 Test Cartridge within 30 minutes of collecting the nasal sample.
T	F	9. This assay should not be used within 30 minutes of administering nasal or throat sprays.
T	F	10. The RNase P internal control has been designed to control for presence of human cellular material in the sample and proper assay execution.
Passing Score is 100% - If you do not pass, reread policy and repeat exam.		Passed
		Failed

<u>Employee</u>	<u>Observer</u>
Printed Name:	Printed Name:
Signature:	Signature: