Intended Use

Sensor and reagent cartridge for use with the IL GEM Premier 5000 analyzer.

Storage Requirement

STORE at 15°-25° C

DO NOT FREEZE

Do not open cartridge protective wrapper until ready to use.

Instructions for Use

To insert a GEM Premier 5000 PAK cartridge:

1. Unpack the cartridge from the protective wrapper. Ensure that the cartridge is dry and within expiration date prior to use.
2. Remove plastic backing from cartridge by pulling tab.
3. Press the Open Door button. The door will release. Fully open the door to the left.
Note: The cartridge must be at 15-25°C prior to use.
4. Position the cartridge with the gray sampling area facing forward. Push the cartridge in until you feel resistance.

Note: The cartridge cannot be inserted all the way into the compartment by hand. A small lip of the cartridge will rest on the door.

5. Guide the system door to the right to close it and move the cartridge into the final position. The door will lock and the warm-up cycle will take approximately 60 minutes. Samples cannot be analyzed during warm-up.

Note: Do not force door closed. If resistance is met when closing the door, please remove the cartridge and restart the process.

6. During the warm-up, the GEM Premier 5000 analyzer will perform Automatic PAK Validation (APV) on all analytes.
7. Run GEM CVP 5 tBili (calibration Valuation Product) to complete the tBili validation process on the GEM Premier 5000 analyzer prior to analyzing patient samples for total Bilirubin.

Refer to the GEM Premier 5000 Operator's Manual for complete instructions.

Warning

WARNING: A bias on certain analytes with capillary samples was observed with RAM Scientific Capillary Tubes (p/n 06 0186) and Fleas (p/n 07 9503). Therefore, do not use RAM Scientific Capillary Tubes (p/n 06 0186) and Fleas (p/n 07 9503) with the GEM Premier 5000. Facilities should evaluate their collection devices prior to clinical use.

GEM 5000 A Solution:

Hazard class: None

Hazard statements: None

Precautionary statements: None

Supplemental hazard information: Up to 2.92% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

GEM 5000 B Solution, GEM 5000 D Solution, GEM 5000 E Solution:

Hazard class: Aquatic chronic 3, H412

Hazard statements: H412: Harmful to aquatic life with long lasting effects.

Precautionary statements: P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: Up to 2.7% of the mixture consists of component of unknown acute toxicity (inhalation) for the human health.

GEM 5000 C Solution:

Warning

Hazard class: Sensitization - Skin, Cat. 1.

Hazard statements: May cause an allergic skin reaction.

Precautionary statements: Avoid breathing vapors/spray. Wear protective gloves. IF ON SKIN: Wash with plenty of water. If skin irritation or a rash occurs: Get medical advice/attention.

Supplemental hazard information: EUH208: Contains m-Phenylendiamine. May produce an allergic reaction. Up to 2.08% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

GEM 5000 Reference Solution:

Warning

Hazard class: Eye Irrit. 2, H319; Aquatic chronic 3, H412.

Hazard statements: H319: Causes serious eye irritation. H412: Harmful to aquatic life with long lasting effects.

Precautionary statements: P280: Wear eye protection/face protection. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313: If eye irritation persists: Get medical advice/ attention. P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: 2.22% of the mixture consists of ingredients of unknown acute toxicity (inhalation) for human health.

GEM 5000 Lysing Solution:

Danger

Hazard class: Eye dam.1, H318; Skin Sens. 1, H317, Aquatic chronic 3, H412.

Hazard statements: H318: Causes serious eye damage H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects.

Precautionary statements: P280: Wear protective gloves/protective clothing/ eye protection/face protection. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P302+P352: IF ON SKIN: Wash with plenty of water. P310: Immediately call a POISON/CENTER/doctor. P273: Avoid release to the environment. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: Contains: Polyethylene glycol octylphenyl ether and Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1). Up to 3.5% of the mixture consists of component of unknown acute toxicity (dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

For *in vitro* diagnostic use.

The GEM® Premier 5000 PAK is protected by one or more of the following US Patents: 5,132,345, 5,286,364, 5,540,828, 6,652,720, 6,872,297, 6,960,466, 7,022,219

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