

MULTISTIX URINE CHEMISTRIES 4, 7 & 10SG on CLINITEK ANALYZER Competency Exam

This exam is to be performed upon initial orientation and training, and semiannual and annual competency testing. (Newly trained staff must perform this test *twice* within the first year.)

**THIS TEST MUST BE PERFORMED IN FRONT OF ANOTHER CERTIFIED OPERATOR
(CERTIFIED OPERATOR = A STAFF MEMBER CURRENTLY APPROVED TO PERFORM TEST)**

****** I HAVE READ THE URINE MULTISTIX ON CLINITEK ANALYZER POLICY AND PROCEDURE BEFORE TAKING
THIS EXAM AS INDICATED BY MY INITIALS HERE: _____.**

WRITTEN TEST QUESTIONS – CIRCLE TRUE OR FALSE	ANSWER	
1. Making up a medical record number (MRN) is okay if it is an emergency.	T	F
2. If I mistrust a result, for any reason, I should repeat the test immediately.	T	F
3. QC on the Clinitek Analyzer is called Siemens Chek-Stix UA Diptube.	T	F
4. QC must be performed daily on the analyzer before performing patient testing.	T	F
5. Only freshly voided (< 1 hour old) urine is an acceptable patient specimen.	T	F
6. Results are shown on analyzer screen for up to 10 minutes.	T	F
7. False positives and false negatives are possible results.	T	F
8. It is important to fully seat the lid onto the Multistix vial, when closing it.	T	F
9. You do not need to wipe away excess urine on the test bar between patient tests and after running QC; contamination can't happen between samples on the analyzer.	T	F
10. If testing manually, all results should be read before 2 minutes.	T	F
You must score 100% on the written test (above) to proceed to test demonstration (below).		SCORE

Printed Name:		License # (Required):		Location:	
Run the Positive Level 2 QC on the analyzer.			EXPIRATION DATE:		
QC LOT #:					
If using Multistix 7 or 10SG, please only list results from the following analytes for this exam:	Leukocytes	Nitrite	Protein	Glucose	

I Performed This Test:		I Observed This Test Being Performed:	
Signature:	Date:	Printed Name:	Date:

Upon Completion: Scan and email a copy to Francis.Corteza@sfdph.org or Cassiusmicho.Santiag@sfdph.org. You may also fax completed tests to POCT Services @ 415.206.3451 or 415.206.3045.