Zuckerberg San Francisco General Hospital POC Quality Control Log CoaguChek® XS Plus (INR)

Liquid Quality Control must be performed:

- 1. With every New lot or New Shipping of Test Strips,
- 2. With every New or New Shipping lot of Quality Control,
- 3. When temperature is out of range, and
- 4. After troubleshooting or preventive maintenance is performed on the meter.

	Lot Number	Exp. Date	Range
Test Strip			N/A
Control 1			
Control 2			

Date	Time	Meter Name	Operator ID	Control 1 Result	Control 2 Result	Corrective Action if Unacceptable/Comments

Reviewed	by:		