

**Zuckerberg San Francisco General Hospital  
Point of Care Testing Temperature Record  
For Room & Refrigerator Temperature Use**

Month:	Year:	Location
--------	-------	----------

Day	Room Temp (*16-25 °C)	Refrigerator Temp (*2-8 °C)	Initials	Comments
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				

**Note: Glucose supplies temperature range is \*16 – 30°C.**

**Reviewed by:** \_\_\_\_\_

\*NOTE: When temperature is outside of defined tolerance limits, corrective actions must be taken and may include, but is not limited to, taking another reading after two hours, submitting a work order to Facilities, and moving contents to another acceptable storage location while problem is being addressed. **Always notify unit supervisor/POCT and document all corrective actions taken.**