

### Cue™ COVID-19 Ag Quality Control Log

Cue™ QC testing is performed:

- When a new shipment or new lot of kits is received
- When problems occur during testing

NOTE: The RNase P internal control has been designed to control for presence of human cellular material in the sample and proper assay execution including sample inhibition, amplification, and assay reagent function. If RNase P is not detected, the Cue COVID-19 Test will return an "Invalid" result.

DATE	CARTRIDGE LOT/EXP	POSITIVE CONTROL LOT/EXP	NEGATIVE CONTROL LOT/EXP	POSITIVE CONTROL	NEGATIVE CONTROL	PASS /FAIL	TESTER'S INITIALS	COMMENTS
EXPECTED RESULTS:				Positive	Negative	PASS		(e.g., new shipment received date)

REVIEWED BY: \_\_\_\_\_ DATE: \_\_\_\_\_